



## Protocol for composing human skin models

### ADMATRIX™

ADMATRIX™ enables you to prepare human skin models easily and handily.

ADMATRIX™ is a highly bioactive ADM (Acellular Dermal Matrix) made from tissue of normal porcine skin, which has the most similar structure to that of human skin. (Fig. 1, 2) Bearing the basement membrane, a higher attachment rate of epithelial cells is shown with ADMATRIX™.

#### ■ How to make human skin models ■

3-Dimensional cultured human skin equivalent with multi-layered epidermis can be composed with ADMATRIX™ used as its scaffold. By using the both side of ADMATRIX™ for culturing, multiple kinds of cells can be incorporated.


#### <Sample procedure to make cultured composite skin: Fig 3>

1. Place ADMATRIX™ set in a culture insert into a 6-well plate, with the dermal side up, and its backings – the filter paper sheets – attached on the both sides. Seed human fibroblasts ( $1 \times 10^5$  cells/mL) onto the surface of ADMATRIX™. SFM medium added with 5% FCS is to be used as culture medium.
2. Incubate for more than one day, and then turn the ADMATRIX™ back over inside the culture insert placed in the 6-well plate. Seed human keratinocytes ( $1 \times 10^5$  cells /mL) onto the basement membrane side of ADMATRIX™, and use SFM medium added with 1% FCS as culture medium.
3. Incubate for 3 days with this medium by liquid phase culture. Then, prepare the culture medium containing SFM and D-MEM at the proportion of 1:1 added with 5% FCS, and with which go on to liquid-phase culture for 24 hours to induce the differentiation of keratinocytes.
4. Exposure the surface of ADMATRIX™ to the air, and continue incubation along with the same culture medium and air-liquid interface culture for 1 week to induce keratinocytes to become multi-layered.




human skin models composed with ADMATRIX™

Fig.1: ADMATRIX™ made from porcine skin tissue, HE stained (x100)



a: Normal porcine skin tissue



b: ADMATRIX™

Fig.2 : ADMATRIX™ porcine type IV collagen, stained (x200)






Fig.3 : human skin model composed from ADMATRIX™ , HE stained (x200)



a: Human keratinocytes (ADMATRIX™, basement membrane side)



b: Human fibroblasts (ADMATRIX™ dermal side)

This product is for RESEARCH USE ONLY. DO NOT use for any other purpose.

Distributor



COSMO BIO CO., LTD.  
Inspiration for Life Science

TOYO EKIMAE BLDG. 2-20, TOYO 2CHOME  
KOTO-KU, TOKYO 135-0016, JAPAN  
TEL : +81-3-5632-9617  
FAX : +81-3-5632-9618  
URL : <http://www.cosmobio.co.jp/>  
e-mail : [export@cosmobio.co.jp](mailto:export@cosmobio.co.jp)

Manufacturer



BCS, Inc

5th floor, Narita Bldg, 3-43-16 Hongo,  
Bunkyo-ku, TOKYO 113-0033 JAPAN  
TEL : +81-3-5840-6531  
FAX : +81-3-5802-5021  
URL : <http://www.bcsinc.co.jp/en/index.htm>  
e-mail : [info@bcsinc.co.jp](mailto:info@bcsinc.co.jp)