

Detection of immune cell checkpoint and functional markers in the tumor microenvironment by the RNA *in situ* hybridization RNAscope[®] assay

Ming-Xiao He, Na Li, Courtney Anderson, Yuling Luo, Nan Su, Xiao-Jun Ma, Emily Park

The RNAscope® assay is a unique RNA ISH technology that identifies RNA expression at the single cell level with morphological context. Here we present the use of RNAscope® for the detection of RNA targets in the tumor microenvironment (TME) that are involved in tumor immunology and immunotherapy:

- · Checkpoint markers
- · Immune cell markers
- · Cytokines and chemokines

Detection of these RNA targets with the RNAscope[®] assay can aid in:

- Localization of specific immune cell types (i.e., cytotoxic lymphocytes and regulatory T cells) in the TME
- Determining spatial relationships between different cell types in the TME
- Characterization of secreted proteins (i.e., cytokines and chemokines)
- Evaluation of immune function in TME beyond enumeration of tumor infiltrating lymphocytes

Introduction

The field of cancer immunotherapy has expanded rapidly in recent years with immune checkpoint inhibitors and other therapeutic approaches such as cancer vaccines and chimeric antigen receptor (CAR) therapy showing promising clinical results. Despite the dramatic and durable responses seen in many patients, our understanding of the immune response to cancer is still limited, and we cannot reliably predict who will or will not benefit from these new interventions. To better stratify patients for immunotherapy treatments, the series of events and biomarkers involved in the cancer-immunity cycle need to be better understood⁽¹⁻³⁾. In addition, spatially mapped expression data at the single-cell level is crucial to understanding the cellular organization and cell-to-cell interactions in the tumor and its complex microenvironment (TME).

RNAscope® is a unique RNA ISH technology that provides single-cell gene expression resolution with spatial and morphological context. The RNAscope® assay detects mRNA and long noncoding RNAs in fresh frozen, fresh fixed, and formalin-fixed paraffin-embedded (FFPE) cells and tissues by utilizing a unique double Z probe design and signal amplification strategy that allows for visualization of target RNA as a single dot, where each dot is an individual RNA molecule⁽⁴⁾. The RNAscope® strategy offers the key benefits of high sensitivity due to the signal amplification method, and high specificity because of the double Z probe design, resulting in a high signal-to-noise ratio in many tissues. In this report we examined 50 selected immune checkpoint and functional markers (summarized in Figure 1) to demonstrate the utility of the RNAscope assay for tumor immunology applications. With the RNAscope® 2.5 HD duplex assay we demonstrate localization of PD-L1 and several immune markers in ovarian and lung cancers (Figures 2-4). Detection of 12 immune checkpoint markers, 14 immune cell markers, and 24 immune function markers, such as cytokines and chemokines, in multiple human tumors is also demonstrated by the RNAscope® brown assay (Figures 5-7; Appendix). The complete data set can be viewed in the appendix available online at

www.acdbio.com/immunotherapy.

Checkpoint	Immune	Cytokines and	
markers	cell markers	chemokines	
BTLA (CD272) CTLA4 IDO1 TIM3 (HAVCR2) LAG3 MICA MICB PD-1 (PDCD1) PD-L1 (CD274) PD-L2 (PDCD1LG2) 4-1BB/CD137 (TNFRSF9) B7-H4 (VTCN1)	CD3 CD8a (CD8A) CD27 CD40 CD40L (CD40LG) CD70 FOXP3 Granulysin (GNLY) ICAM1 ICOS Perforin (PRF1) OX40L (TNFSF4) OX40 (TNFRSF4) HVEM (TNFRSF14)	CCL5 CCL7 CX3CL1 CX3CR1 CXCL9 CXCL10 CXCL13 IFNB1 IFNB1 IFNY (IFNG) IL1B IL4 IL6 IL7	IL10 IL12A IL13 IL17A IL18 IL20 IL21 IL33 TGFβ (TGFB1) TNFα (TNFA) VEGF-A (VEGFA)

FIGURE 1. Immune markers examined in this study by the RNAscope® assay. Fifty immune system markers were examined in this study, including 12 immune checkpoint markers, 14 immune cell markers, and 24 cytokines and chemokines.



FIGURE 2. Simultaneous detection of *PD-L1* and *CD8a* mRNAs. The manual RNAscope[®] 2.5 HD Duplex assay was used to detect expression of the immune checkpoint marker *PD-L1* (green chromogen) and the immune cell marker *CD8a* (red chromogen) in human ovarian (A) and lung (B) tumors. Inset shows enlarged region outlined in white. Asterisk denotes *PD-L1+/CD8-* cells; arrowhead denotes *PD-L1+/CD8+* cells. Note that in the lung cancer sample (B) *PD-L1* is profusely expressed in both tumor cells and immune cells, with abundant *CD8+* immune cell infiltration. However, in the ovarian cancer sample (A) *PD-L1* expression is primarily restricted to the immune cells. S, stroma; T, tumor. White line delineates tumor margin. 40x magnification.

Ovarian cancer

Lung cancer



FIGURE 3. Simultaneous detection of *IFN*_Y and immune cell markers *CD3*, *CD8a* and *CD4* mRNAs. The manual RNAscope[®] 2.5 HD Duplex assay was used to detect expression of the immune function marker *IFN*_Y and the immune cell markers *CD3* (A-B), *CD8a* (C-D), or *CD4* (E-F) in human ovarian (A, C, E) and lung (B, D, F) tumors. *CD3*, *CD8a*, and *CD4* were detected using green chromogen and *IFN*_Y was detected using red chromogen. Inset shows enlarged region outlined in white. Asterisk denotes *CD+/IFN*_Y- cells; arrowhead denotes *CD+/IFN*_Y+ cells; arrow denotes *CD-/IFN*_Y+ cells. S, stroma; T, tumor. White line delineates tumor margin. 40x magnification.



FIGURE 4. Simultaneous detection of *FOXP3* and *CD4* mRNAs. The manual RNAscope[®] 2.5 HD Duplex assay was used to detect expression of the immune cell markers *CD4* (green chromogen) and *FOXP3* (red chromogen) in human ovarian (A) and lung (B) tumors. Inset shows enlarged region outlined in white. Asterisk denotes *CD4+/FOXP3*- cells; arrowhead denotes *CD4+/FOXP3*+ regulatory T cells. S, stroma; T, tumor. White line delineates tumor margin. 40x magnification.



FIGURE 5. Detection of immune checkpoint markers in human tumors. The RNAscope® manual brown assay was used to detect expression of several immune checkpoint markers in multiple human tumors. 40x magnification.



FIGURE 6. Detection of immune cell markers in human tumors. The RNAscope® manual brown assay was used to detect expression of several immune cell markers in multiple human tumors. 40x magnification.



FIGURE 7. Detection of cytokines, chemokines, and their receptors in human tumors. The RNAscope® manual brown assay was used to detect expression of several cytokines, chemokines, and their receptors in multiple human tumors. 40x magnification.

Conclusions

Detecting RNA biomarker expression at the single-cell level while preserving spatial information is critical to understanding cellular organization and cell-to-cell interactions in the cancer-immunity cycle. Here we show that the RNAscope® assay is able to detect immune checkpoint markers and immune function markers in a variety of human tumors. Overall these results demonstrate

Markers	Probe name	Cat No.		
Checkpoint markers				
B7-H4 (VTCN1)	Hs-VTCN1	418081		
BTLA (CD272)	Hs-BTLA	401601		
CD137/4-1BB (TNFRSF9)	Hs-TNFRSF9	415171		
CTLA4	Hs-CTLA4	554341		
ID01	Hs-ID01	602681		
TIM3 (HAVCR2)	Hs-HAVCR2	560681		
LAG3	Hs-LAG3	553931		
MICA	Hs-MICA	427161		
MICB	Hs-MICB	427181		
PD-1 (PDCD1)	Hs-PDCD1	602021		
PD-L1 (CD274)	Hs-CD274	600861		
PD-L2 (PDCD1LG2)	Hs-PDCD1LG2	551891		
Immune cell markers				
CD3	Hs-CD3-pool	426621		
CD8a (CD8A)	Hs-CD8A	560391		
CD27	Hs-CD27	415451		
CD40	Hs-CD40	578471		
CD40L (CD40LG)	Hs-CD40LG	542341		
CD70	Hs-CD70	419331		
FOXP3	Hs-FOXP3	418471		
Granulysin (GNLY)	Hs-GNLY	407371		
HVEM (TNFRSF14)	Hs-TNFRSF14	319731		
ICAM1	Hs-ICAM1	402951		
ICOS	Hs-ICOS	407141		
OX40L (TNFSF4)	Hs-TNFSF4	427201		
OX40 (TNFRSF4)	Hs-TNFRSF4	412381		
Perforin (PRF1)	Hs-PRF1	407381		

the utility of the RNAscope[®] technology in studying tumor immunology and key targets of immunotherapy in the tumor and its complex microenvironment.

Appendix containing all data discussed in this study is available at acdbio.com/immunotherapy

Markers	Probe name	Cat No.		
Cytokines and chemokines				
CCL5	Hs-CCL5	549171		
CCL7	Hs-CCL7	425261		
CX3CL1	Hs-CX3CL1	411261		
CX3CR1	Hs-CX3CR1	411251		
CXCL9	Hs-CXCL9	440161		
CXCL10	Hs-CXCL10	311851		
CXCL13	Hs-CXCL13	311321		
IFNB1	Hs-IFNB1	417071		
IFN-γ (IFNG)	Hs-IFNG	310501		
IL1B	Hs-IL1B	310361		
IL4	Hs-IL4	315191		
IL6	Hs-IL6	310371		
IL7	Hs-IL7	424251		
IL10	Hs-IL10	602051		
IL12A	Hs-IL12A	402061		
IL13	Hs-IL13	586241		
IL17A	Hs-IL17A	310931		
IL18	Hs-IL18	400301		
IL20	Hs-IL20	412201		
IL21	Hs-IL21	401251		
IL33	Hs-IL33	400111		
TGF-β (TGFB1)	Hs-TGFB1	400881		
TNF-α (TNFA)	Hs-TNFA	310421		
VEGF-A (VEGFA)	Hs-VEGFA	423161		

Probes for duplex assay				
Markers	Probe name	Cat #		
PD-L1/CD8A	Hs-CD274	600681		
	Hs-CD8A-C2	560391-C2		
PD-L1/CD3*	Hs-CD274	600681		
	Hs-CD3-pool-C2	426621-C2		
CD3/IFNG	Hs-CD3-pool	426621		
	Hs-IFNG-C2	310501-C2		
CD8a/IFNG	Hs-CD8A	560391		
	Hs-IFNG-C2	310501-C2		
CD4/IFNG	Hs-CD4	605601		
	Hs-IFNG-C2	310501-C2		
CD4/FOXP3	Hs-CD4	605601		
	Hs-FOXP3-C2	418471-C2		

* Data not included in this report.

References

- 1. Pardoll DM. The blockade of immune checkpoints in cancer immunotherapy. Nat Rev Cancer. 2012; 12(4): 252-264.
- Chen DS and Mellman I. Oncology meets immunology: the cancer-immunity cycle. Immunity. 2013; 39(1): 1-10.
- Tumeh PC et al. PD-1 blockade induces new responses by inhibiting adaptive immune response. Nature. 2014; 515(7528): 568-571.
- Wang F, Flanagan J, Su N, Wang LC, Bui S, Nielson A, Wu X, Vo HT, Ma XJ, Luo Y. RNAscope: A novel in situ RNA analysis platform for formalin-fixed, paraffin-embedded tissues. J Mol Diagn. 2012; 14(1):22-9.

DISCLAIMER

All samples used in this study have been qualified for RNA integrity with *PPIB* and *dapB* control probes before target probe analysis. Always qualify your samples as describe in the getting started section www.acdbio.com/technical-support/getting-started. You may observe a difference in staining pattern, due to variation in tissue preparation methods or associated biology within the tissue samples.

To read additional Application Notes, visit **www.acdbio.com**

For Research Use Only. Not for diagnostic use. RNAscope is a registered trademark of Advanced Cell Diagnostics, Inc. in the United States or other countries. All rights reserved. ©2015 Advanced Cell Diagnostics, Inc. Doc #: 51-056/RevB/Effective Date03/31/2016

