

## Contact Information

70-262 DRx\_UM\_PHeX\_TOTAL\_0914V5



**DiscoverRx Corporation**  
(World Wide Headquarters)  
42501 Albrae Street  
Fremont, CA 94538  
United States

t | 1.510.979.1415  
f | 1.510.979.1650  
toll-free | 1.866.448.4864

**DiscoverRx Corporation Ltd.**  
(Europe Headquarters)  
Faraday Wharf, Holt Street  
Birmingham Science Park Aston  
Birmingham, B7 4BB  
United Kingdom

t | +44.121.260.6142  
f | +44.121.260.6143

**KINOMEScan®**  
*A division of DiscoverRx*  
11180 Roselle Street, Suite D  
San Diego, CA 92121  
United States

t | 1.800.644.5687  
f | 1.858.630.4600

**BioSeek®**  
*A division of DiscoverRx*  
310 Utah Avenue, Suite 100  
South San Francisco, CA 94080  
United States

t | 1.650.416.7600  
f | 1.650.416.7625

[www.discoverx.com](http://www.discoverx.com)



# PathHunter® eXpress Total GPCR Internalization Assays

## User Manual

**CONTENTS**

LEGAL SECTION	PAGE 3
INTENDED USE	PAGE 4
TECHNOLOGY PRINCIPLE	PAGE 4
PROTOCOL OVERVIEW	PAGE 5
KIT CONTENTS AND STORAGE CONDITIONS	PAGE 6
MATERIALS PROVIDED	PAGE 7
ADDITIONAL MATERIALS REQUIRED (NOT PROVIDED)	PAGE 7
RECOMMENDED MATERIALS	PAGE 7
ASSAY INCUBATION AND CELL PLATING REAGENT REQUIREMENTS	PAGE 8
COMPOUND PREPARATION AND INCUBATION TEMPERATURES	PAGE 8
USE OF PLASMA OR SERUM CONTAINING SAMPLES	PAGE 8
THAWING AND PLATING FROZEN CELLS	PAGE 9
ASSAY PROCEDURE — AGONIST DOSE RESPONSE	
◆ PLATE MAP	PAGE 10
◆ PROTOCOL	PAGE 10
◆ QUICK START PROCEDURE	PAGE 12
ASSAY PROCEDURE — ANTAGONIST DOSE RESPONSE	
◆ PLATE MAP	PAGE 13
◆ PROTOCOL	PAGE 13
◆ QUICK START PROCEDURE	PAGE 16
FREQUENTLY ASKED QUESTIONS	PAGE 17
APPENDIX A: RELATED PRODUCTS	PAGE 20

**NOTES:**

## NOTES:

## LEGAL SECTION

**This product and/or its use is covered by one or more U.S. patents #7,135,325 B2, #8,101,373 B2, and/or foreign patent applications and trade secrets that are either owned by or licensed to DiscoverX® Corporation.**

### LIMITED USE LICENSE AGREEMENT

The designated cells and reagents purchased from DiscoverX are restricted in their use. DiscoverX has developed an assay for translocation and internalization ("Assay") employing genetically modified cells ("Cells") and detection reagents ("Reagents") (collectively referred to as "Materials"). The Cells and Reagents are designed and optimized to be used together in the Assay. DiscoverX wishes to ensure that these Cells and Reagents are used properly and effectively. By purchasing the Materials you recognize and agree to the restrictions.

1. The Materials are not transferable and will be used only at the site for which they were purchased. Transfer to another site owned by Purchaser will be permitted only upon written request by Purchaser followed by subsequent written approval by DiscoverX.
2. Purchaser will not analyze the Reagents nor have them analyzed on Purchaser's behalf.
3. Purchaser will use only the Reagents supplied by DiscoverX or an authorized DiscoverX distributor for the Assays.

If the purchaser is not willing to accept the limitations of this limited use statement and/or has any further questions regarding the rights conferred with purchase of the Materials, please contact:

**Licensing Department**  
DiscoverX Corporation  
42501 Albrae Street  
Fremont, CA 94538 USA  
tel | 1.510.979.1415 x104  
[info@discoverx.com](mailto:info@discoverx.com)

**For some products/cell lines, certain 3rd party gene specific patents may be required to use the cell line. It is the purchaser's responsibility to determine if such patents or other intellectual property rights are required.**

## INTENDED USE


**PathHunter® eXpress Total GPCR Internalization Assays** are ready to use complete kits that contain everything you need to measure internalized GPCRs in live cells but without imaging, antibodies or radioactivity. The eXpress kits include single use vials of frozen cells stably expressing the GPCR of interest, optimized cell plating reagent, chemiluminescent detection reagents and plates\*. Simply thaw and plate the pre-validated cells and challenge with compound 24 or 48 hours later. Whether you are studying receptor recycling and kinetics, identifying novel inhibitors of receptor internalization, or determining mechanism of action of your lead compounds, the ready-to-assay eXpress format eliminates the need for lengthy, expensive and time consuming cell culture and makes functional testing fast and convenient. Assays are designed for both 96-well plate analyses and kits include enough cells and detection reagents for either 100, 200 or 1,000 data points.

\*Test compounds are not included and must be provided by the researcher.

## NOTES:

## TECHNOLOGY PRINCIPLE

**PathHunter Total GPCR Internalization Assays** provide a direct and quantitative measurement of internalized GPCR protein localized in early endosomes using  $\beta$ -galactosidase ( $\beta$ -gal) enzyme fragment complementation (EFC, Figure 1). These are available in one of two formats A) The small, 42 amino acid enzyme fragment of  $\beta$ -gal called ProLink™ (PK) is fused to the GPCR of interest and the larger, complementing enzyme fragment termed Enzyme Acceptor, or EA, is localized to the endosomes and B) The small, 42 amino acid enzyme fragment of  $\beta$ -gal called Pro-Link (PK) is localized to the endosomes and the larger, complementing enzyme fragment termed Enzyme Acceptor, or EA is fused to the GPCR of interest. GPCR activation results in internalization of the receptor in endosomes. This action forces complementation of the two enzyme fragments, resulting in an increase in enzyme activity that is easily measured using chemiluminescent PathHunter Detection Reagents.



**Figure 1. PathHunter® Total GPCR Internalization Assay Principle.** Activation of the GPCR results in internalization of the receptor in endosomes and formation of a functional  $\beta$ -gal enzyme capable of hydrolyzing substrate and generating chemiluminescent signal.

APPENDIX A: RELATED PRODUCTS


Description	Catalog Number	For more information, visit:
Control Ligands	Many	<a href="http://www.discoverx.com/ligands/control_ligands.php">www.discoverx.com/ligands/control_ligands.php</a>
AssayComplete™ Cell Plating Reagents	93-0563R0A 93-0563R5A 93-0563R7A 93-0563R28A	<a href="http://www.discoverx.com/cell/cell_plating_reagents.php">http://www.discoverx.com/cell/cell_plating_reagents.php</a>
PathHunter® eXpress β-Arrestin GPCR Assays	Many	<a href="http://www.discoverx.com/gpcrs/express_arrestin.php">www.discoverx.com/gpcrs/express_arrestin.php</a>
PathHunter® eXpress β-Arrestin Orphan GPCR Assays	Many	<a href="http://www.discoverx.com/gpcrs/express_orphan.php">www.discoverx.com/gpcrs/express_orphan.php</a>
PathHunter® eXpress β-Arrestin Ortholog GPCR Assays	Many	<a href="http://www.discoverx.com/gpcrs/express_ortholog.php">www.discoverx.com/gpcrs/express_ortholog.php</a>

PROTOCOL OVERVIEW

Please read the entire protocol completely before running the assay. Successful results depend on performing these steps correctly. Refer to the cell-line specific datasheet for additional information on optimized cell plating reagent and reference ligand. For additional information or Technical Support, contact DiscoverX or visit [www.discoverx.com](http://www.discoverx.com).

The following steps are required to monitor the fate of Total and internalized GPCRs using a PathHunter eXpress Total GPCR Internalization assay (Figure 2).

1. Thaw and plate frozen, assay-ready eXpress cells (page 9).
2. Dilute and add compounds.
3. Perform functional assay in agonist (page 10) or antagonist mode (page 13).



**Figure 2.** Monitor functional GPCR responses to compound challenge using the fast and simple PathHunter eXpress Total GPCR Internalization assay protocol.  
\*Please refer to datasheet for any variations in assay conditions.

## KIT CONTENTS AND STORAGE CONDITIONS

**PATHHUNTER EXPRESS TOTAL GPCR INTERNALIZATION KIT COMPONENTS REQUIRE MULTIPLE STORAGE TEMPERATURES. OPEN BOXES IMMEDIATELY AND STORE CONTENTS AS INSTRUCTED.**

### **SHELF LIFE:**

Use kit within 6 months from the date of receipt under proper storage conditions.

### **Box 1: PATHHUNTER EXPRESS TOTAL GPCR INTERNALIZATION CELLS:**

#### **STORAGE:**

Short term (2 weeks or less): **Store vials -80°C immediately upon arrival.**

Long term (greater than 2 weeks): **Place vials in the vapor phase of liquid nitrogen (N<sub>2</sub>).**

PathHunter eXpress Total GPCR Internalization cells arrive frozen on dry ice. Cells are delivered in individual vials containing 1x10<sup>6</sup> cells in 100 µL of freezing medium. Each vial contains sufficient cell numbers to generate (1) 96-well microplate prepared at the seeding density described.

When removing cryovials from liquid N<sub>2</sub> storage, use tongs and place immediately on dry ice in a covered container. Wait at least one minute for any liquid N<sub>2</sub> inside the vial to evaporate and proceed with the thawing protocol (page 9). **Do not touch the bottom of the tubes at any time to avoid inadvertent thawing of the cells. If cells are not frozen upon arrival, do not proceed. Contact technical support.**

### **Box 2: PATHHUNTER DETECTION REAGENT AND CP REAGENT: Store at -20°C**

Once thawed, store the Cell Plating (CP) Reagent at 4°C. Avoid multiple freeze/thaw cycles. In rare instances, the CP Reagent may be yellow in color after thawing. Although this indicates a slight change in pH, continue with the assay as this does not impact assay performance.

Thaw the PathHunter Detection Reagents at room temperature before use, and after thawing, store reagents for up to 7 days at 4°C. The reagents can tolerate up to three freeze-thaw cycles with no impact on performance. Once made, the working solution is stable for 24 hours at room temperature.

### **Box 3: 96-WELL TISSUE CULTURE TREATED PLATES: Store at Room Temperature**

## FREQUENTLY ASKED QUESTIONS (CONTINUED)

**Q: What is the length of compound incubation required for optimal detection of internalization and recycling of the receptor?**

A: Optimal compound incubation times are somewhat target specific. However, the majority of receptors plateau within 2 to 3 hours. Therefore, the PathHunter GPCR Internalization eXpress assays were developed using a single, universal protocol that includes a 3 hour compound incubation step.

**Q: Is the EFC signal impacted by a change in pH?**

A: No. The enzyme fragments are always in the cytosol, so the pH does not change.

## FREQUENTLY ASKED QUESTIONS (CONTINUED)

**Q: Why do longer incubation times with Detection Reagents lead to a higher signal?**

A: The complemented  $\beta$ -galactosidase ( $\beta$ -gal) enzyme is continually turning over the substrate over time. Theoretically, the signal continues to increase until the substrate is exhausted. Therefore, the longer you incubate the reaction, the higher the RLU values.

**Q: What is the shelf life of the eXpress kits?**

A: We recommend that eXpress kits should be used within 6 months of receipt under proper storage conditions. For short term (2 weeks or less), store eXpress cells at  $-80^{\circ}\text{C}$ . For long term storage (more than 2 weeks), store in the vapor phase of liquid nitrogen ( $\text{N}_2$ ). Store the Detection Reagent Kit at  $-20^{\circ}\text{C}$ . Refer to the kit label for lot specific expiration date information.

**Q: What if my Cell Plating Reagent changes from a red/pink color to yellow after freezing/thawing?**

A: If the Cell Plating Reagent changes color from red/pink to yellow after thawing, please continue with the assay according to the product insert. We have observed this color change on rare occasions and have confirmed that it will not affect assay performance.

**Q: Can I use this assay to test human plasma or serum samples?**

A: Yes. PathHunter eXpress GPCR assays tolerate up to 80% serum or plasma. First, prepare a standard curve of spiked ligand in neat, heparinized plasma (or mouse, human serum). Add samples directly to the cells (no further dilution – 100% plasma in the well). After stimulation, remove the plasma or serum sample and replace with fresh CP reagent before addition of the PathHunter Detection Reagents. It has been shown that EDTA anti-coagulated plasma inhibits EFC and should be avoided for these types of studies.

**Q: What instruments can I use to read the plates?**

A: Any bench top luminometer will work with the PathHunter eXpress GPCR Assays. Below is a partial list of commercially available luminometers that have been used to validate our assays:

**Turner Biosystems:** Modulus Microplate

**GE Healthcare Life Sciences:** LEADseeker™, FarCyte™

**BMG Labtech:** PHERAstar Plus, LUMIstar Omega

**Perkin Elmer:** TopCount®, VICTOR II or V, Fusion, LumiCount, EnVision, MicroBeta® (Trilux), ViewLux

**Molecular Devices:** CLIPR™, LJI Acquest, LJI Analyst, LJI Analyst HT, LJI Analyst GT, Gemini, SpectraMax®, Flexstation™, LMax

**Tecan:** Ultra Evolution

**Beckman Coulter** – CRI

**Berthold Technologies:** Mithras LB 940

**Hamamatsu:** FDSS6000, FDSS/RayCatcher

## MATERIALS PROVIDED

Description	Contents			Storage
<b>Box 1:</b> PathHunter eXpress Total GPCR Internalization Cells	1 vial 1x10 <sup>6</sup> cells ea	2 vials 1x10 <sup>6</sup> cells ea	10 vials 1x10 <sup>6</sup> cells ea	$-80^{\circ}\text{C}$ (short) Liquid $\text{N}_2$ (long)
<b>Box 2:</b> PathHunter Detection Reagents - Cell Assay Buffer - Substrate Reagent 1 - Substrate Reagent 2* Cell Plating Reagent <sup>†</sup>	100dp  5.7 mL 1.5 mL 0.3 mL 1 X 20.0 mL	200 dp  9.5 mL 2.5 mL 0.5 mL 2 X 20.0 mL	1,000 dp  57.0 mL 15.0 mL 3.0 mL 2 X 100 mL	$-20^{\circ}\text{C}$
<b>Box 3:</b> 96-well Tissue Culture Treated Plates	1 plate	2 plates	10 plates	Room Temp

<sup>†</sup>Refer to cell-line specific data sheets for optimized Cell Plating Reagent included with each kit.

\*Centrifuge vial before opening to maximize recovery.

## ADDITIONAL MATERIALS REQUIRED (NOT PROVIDED)

The following additional materials are required but not provided:

1. Pipettes and pipette tips
2. Tissue culture disposables
3. GPCR control agonist as recommended in the cell line specific datasheet. Visit [www.discoverx.com/pathway\\_assays/control\\_ligands.php](http://www.discoverx.com/pathway_assays/control_ligands.php) for the complete DiscoverX offering.
4. GPCR Test compound(s)
5. Disposable Reagent Reservoir such as Thermo Scientific, Cat. #8094 or similar
6. 96-well V-bottom compound dilution plates (DiscoverX, Cat. #92-0011)
7. Multi-mode or luminescence plate reader.

## RECOMMENDED MATERIALS

The following products\* are recommended:

- CytoTracker™ LDH Quantification Kit (DiscoverX, Cat. # 92-2002)
- CytoTracker™ Glutathione Quantification Kit (DiscoverX, Cat. # 92-2003)
- CytoTracker™ DNA Damage Quantification Kit (DiscoverX, Cat. # 92-2004M)

\* Products not available in all countries. Please inquire.

## ASSAY INCUBATION AND CELL PLATING REAGENT REQUIREMENTS

Each PathHunter eXpress Total GPCR Internalization Assay has been validated for optimal assay performance at either 24 or 48 hours post-thaw. Although most targets perform similarly at both time points, for optimal assay performance we recommend you perform the assay according to the protocol provided in the cell line specific datasheet using both the recommended time point and CP Reagent. **Always use the CP Reagent included in the kit and DO NOT substitute from an alternate kit at any time.**

### NOTE:

Use special caution when testing multiple targets in the same experiment as targets may have different incubation times and CP Reagent requirements.

## COMPOUND PREPARATION AND INCUBATION TEMPERATURE

PathHunter eXpress Total GPCR Internalization Assays are routinely carried out in the presence of  $\leq 1\%$  solvent (i.e. DMSO, ethanol, PBS or other). As solvents can affect assay performance, optimize the assay conditions accordingly if other solvents or solvent concentrations are required.

To validate each PathHunter eXpress Total GPCR Internalization Assay, reference ligand was diluted in the recommended CP Reagent containing appropriate solvent. For preparation of test compounds, we recommend preparing the dilutions using the CP Reagent provided in the kit. For antibodies or other compounds that may be sensitive to serum and/or other assay components, dilutions can be prepared in either Hanks Buffered Salt Solution (HBSS) + 10 mM HEPES 0.1% Bovine Serum Albumin (BSA) or OptiMEM® + 0.1% (BSA) without affecting assay performance.

The kinetics of ligand-induced receptor internalization can vary depending on the target and temperature used during the compound incubation step. For optimal assay performance, we recommend you perform the compound incubation step according to the protocol provided in the cell line specific datasheet. **Always use the incubation temperature recommended for the kit you are testing.**

## USE OF PLASMA OR SERUM CONTAINING SAMPLES

PathHunter eXpress Total GPCR Internalization Assay can be run in the presence of high levels of serum or plasma without negatively impacting assay performance. Standard curves of control ligand can be prepared in neat, heparinized plasma and added directly to the cells (without further dilution, i.e. 100% plasma in the well). After ligand stimulation, the samples should be removed and replaced with fresh CP Reagent before the addition of the PathHunter Detection Reagents.

### NOTE:

EDTA anti-coagulated plasma samples do not give a positive response in the assay. Therefore, the choice of anti-coagulant treatment is very important.

## FREQUENTLY ASKED QUESTIONS

### Q: I did not see a signal with my control agonist.

A: There may be differences in agonist purchased from different vendors. Confirm that the control agonist used is the same ligand used in the dose response shown in the provided cell-specific data sheet.

### Q: Can the source of my agonist or antagonist compound impact my assay performance?

A: Yes, the vendor/source of compound can impact assay performance dramatically. Compounds can vary in purity from vendor to vendor. In addition, vendors will recommend different diluents (methanol, NaOH, ethanol, DMSO, water), different treatments (boiling, freeze/thaw, etc) or different storage temperatures for the same compound. Each PathHunter eXpress target has been QC tested and validated using a reference ligand. Information on the reference ligand used for each assay (including the vendor source and catalog number) can be found on the cell line specific datasheet. For optimal assay performance, we recommend using control ligands provided by DiscoverX. Visit [www.discoverx.com/ligands/control\\_ligands](http://www.discoverx.com/ligands/control_ligands) for the complete DiscoverX offering.

### Q: I did not see a response with my compound.

A1: The concentration of DMSO or Ethanol used for dilution is too high. Maintain concentration of the agonist/antagonist diluent at  $\leq 1\%$ .  
A2: Confirm that the final ligand concentration is correct. Some ligands are "sticky" and difficult to dissolve.  
A3: Confirm that the cell line responds to the control agonist.  
A4: Repeat the experiment using a new lot of control agonist.

### Q: My cells arrived thawed. Can I use them?

A: No. Call technical support for a replacement.

### Q: How long is the prepared detection reagent good for?

A: The working detection reagent solution must be used within 8 hours of mixing.

### Q: How long is the signal stable for?

A: The signal is stable for 24 hours after addition of detection reagent.

### Q: My cells are floating after the 48 hours incubation.

A: The cells are not viable, contact technical support for a replacement.

### Q: Can I switch plates or should I use the plate provided?


A: You can use any clear bottom white or opaque walled plate.

### Q: What if cells are not completely adherent after 24/48 hrs incubation?

A: For certain targets, cells may not be completely adherent after 24 hours, but still greater than 80% viable. Please continue on with the protocol as described in the product insert.



## QUICK-START PROCEDURE: ANTAGONIST DOSE RESPONSE



\*Please refer to the cell line specific datasheet for any variation in assay conditions.


## THAWING AND PLATING FROZEN CELLS

The following steps outline the procedure for thawing and plating frozen PathHunter eXpress Total GPCR Internalization cells from freezer vials:

1. Pre-warm CP Reagent in a 37°C water bath.
2. Remove cell vial(s) from -80°C or liquid N<sub>2</sub> vapor phase storage and place immediately on dry ice prior to thawing. **DO NOT EXPOSE VIALS TO ROOM TEMPERATURE.**  
**NOTE:**  
When removing cryovials from liquid N<sub>2</sub>, place immediately on dry ice in a covered container. Wait at least one minute before opening for any liquid N<sub>2</sub> inside the vial to evaporate.
3. Place the cell vial(s) **briefly** (10 seconds to 1 min) in a 37°C water bath until only small ice crystals remain and the cell pellet(s) is almost completely thawed.
4. Add 0.5 mL of pre-warmed CP Reagent to the cell vial. Pipette up and down gently several times to ensure that the cells are evenly distributed.
5. Immediately transfer the cells to 11.5 mL of pre-warmed CP Reagent and pour into a disposable reagent reservoir.
6. Plate 100 µL of cells into each well of the 96-well tissue culture plate.
7. After seeding the cells into the microplate, incubate for either 24 or 48 hours at 37°C, 5% CO<sub>2</sub>. Please refer to the cell line specific datasheet for the incubation time required for the PathHunter eXpress Total GPCR Internalization Assay you are testing.

## ASSAY PROCEDURE - AGONIST DOSE RESPONSE

The steps outlined below provide the assay volumes and procedure for performing agonist assays using the PathHunter eXpress Total GPCR Internalization cells and PathHunter Detection Reagents. Although plate layouts and experimental designs may vary, we recommend performing a 12-point dose curve for each compound using at least duplicate wells for each dilution. The protocol and volumes described below are designed for a complete 96-well plate.



## DAY 2 OR 3: AGONIST COMPOUND PREPARATION AND ADDITION

1. Dissolve agonist compound in the vehicle of choice (DMSO, Ethanol, PBS or other) at the desired concentration.
2. Prepare 3-fold serial dilutions of agonist compound in CP Reagent containing the appropriate solvent (DMSO, ethanol, PBS or other). The concentration of each dilution should be prepared at **11X** of the final screening concentration (i.e. 10  $\mu$ L compound + 100  $\mu$ L of cells). For each dilution, the final concentration of solvent should remain constant.

### Preparation of 12-point dose curve serial dilutions:

We recommend starting with a concentration that is **50X** the expected  $EC_{50}$  value for the compound (**550X**  $EC_{50}$  would be the final working concentration).

**Example:** If the expected  $EC_{50}$  is 10 nM, prepare the highest working concentration at 5.5  $\mu$ M.

- a. For each compound tested, label tubes 1 through 12.
- b. Add 60  $\mu$ L of CP Reagent containing appropriate solvent to tubes #1-11.
- c. Prepare a working concentration of agonist compound in appropriate CP Reagent.
- d. Add 90  $\mu$ L of the working concentration of agonist compound to tube #12.
- e. Remove 30  $\mu$ L of diluted compound from tube #12, add it to tube #11 and mix gently by pipetting up and down. Discard the pipet tip.

## SUBSTRATE PREPARATION AND ADDITION

1. During the incubation period, prepare a working stock of PathHunter Detection Reagents by mixing **19 parts** Cell Assay Buffer, **5 parts** Substrate Reagent 1 and **1 part** Substrate Reagent 2.

Component	Entire Plate (96 wells)
Cell Assay Buffer	4.75 mL
Substrate Reagent 1	1.25 mL
Substrate Reagent 2	0.25 mL

### NOTE:

The working solution is stable for up to 8 hours at room temperature.

2. Add 55  $\mu$ L of prepared detection reagent per well and incubate for 60 minutes at room temperature (23°C). **DO NOT pipette up and down in the well to mix or vortex/shake plates.**
3. Read samples on any standard luminescence plate reader.
4. Use GraphPad Prism® or other comparable program to plot your antagonist dose response.

- a. Label tubes 1 through 12.
  - b. Add 60  $\mu$ L of CP Reagent containing appropriate solvent to tubes #1-11.
  - c. Prepare a working stock of antagonist compound in the appropriate CP Reagent.
  - d. Add 90  $\mu$ L of the working concentration of antagonist compound to tube #12.
  - e. Remove 30  $\mu$ L of diluted compound from tube #12, add it to tube #11 and mix gently by pipetting up and down. Discard the pipet tip.
  - f. With a clean pipet tip, remove 30  $\mu$ L of diluted compound from tube #11, add it to the tube #10 and mix gently by pipetting up and down. Discard the pipet tip.
  - g. Repeat this process 7 more times, preparing serial dilutions from right to left across the plate. **DO NOT add antagonist compound to tubes #1 and 2.** Add only CP Reagent containing appropriate solvent. These samples serve as the no antagonist controls and complete the dose curve.
  - h. Repeat process when testing additional compounds.
  - i. Set compounds aside until antagonist compounds are ready to be added.
3. Remove PathHunter<sup>®</sup> eXpress Total GPCR Internalization cells (previously plated on day 1) from the incubator.
  4. Transfer 5  $\mu$ L from tubes #1-12 to each well according to the plate map on page 13.
  5. Incubate for 30 minutes at 37°C.

#### AGONIST COMPOUND PREPARATION AND ADDITION

1. During the antagonist incubation, determine the EC<sub>80</sub> concentration of the agonist from the agonist dose response curve (described on pages 10-12). Prepare a **22X** EC<sub>80</sub> concentration of agonist compound as shown below:  
**Example:** If the EC<sub>80</sub> of the agonist compound is 10 nM, prepare a stock at 220 nM.
2. Add 5  $\mu$ L of agonist compound to each well. Add 5  $\mu$ L of CP Reagent to the no agonist wells (column 1).
3. Incubate for 3 hours @ 25°C or 37°C\*.

**NOTE:**

\*Please refer to the cell line specific datasheet for any variation in assay conditions.

- f. With a clean pipet tip, remove 30  $\mu$ L of diluted compound from tube #11, add it to tube #10 and mix gently by pipetting up and down. Discard the pipet tip.
  - g. Repeat this process 8 more times, preparing serial dilutions from right to left across the tubes. **DO NOT add agonist compound to tube #1.** Add only appropriate CP Reagent containing appropriate solvent. This sample serves as the no agonist control and completes the dose curve.
  - h. Repeat this process for each compound to be tested.
  - j. Set compounds aside until agonist compounds are ready to be added.
3. Remove PathHunter<sup>®</sup> eXpress Total GPCR Internalization cells (previously plated on day 1) from the incubator.
  4. Transfer 10  $\mu$ L from tubes #1-12 to each well according to the plate map on p.10.
  5. Incubate for 3 hours @ 25°C or 37°C\*.

**\*NOTE:**

Please refer to the cell line specific datasheet for the compound incubation temperature required for the PathHunter<sup>®</sup> eXpress Total GPCR Internalization Assay you are testing.

#### SUBSTRATE PREPARATION AND ADDITION

1. During the incubation period, prepare a working stock of PathHunter Detection Reagents by mixing **19 parts** Cell Assay Buffer, **5 parts** Substrate Reagent 1 and **1 part** Substrate Reagent 2.


Component	Entire Plate (96 wells)
Cell Assay Buffer	4.75 mL
Substrate Reagent 1	1.25 mL
Substrate Reagent 2	0.25 mL

**NOTE:**

The working solution is stable for up to 8 hours at room temperature.

2. Add 55  $\mu$ L of prepared detection reagent per well and incubate for 60 minutes at room temperature (23°C). **DO NOT pipette up and down in the well to mix or vortex/shake plates.**
3. Read samples on any standard luminescence plate reader.
4. Use GraphPad Prism<sup>®</sup> or other comparable program to plot your agonist dose response.


## QUICK-START PROCEDURE: AGONIST DOSE RESPONSE



\*Please refer to the cell line specific datasheet for any variation in assay conditions.

## ASSAY PROCEDURE — ANTAGONIST DOSE RESPONSE

The steps outlined below provide the assay volumes and procedure for performing antagonist assays using the PathHunter eXpress Total GPCR Internalization cells and PathHunter Detection Reagents. Although plate layouts and experimental designs may vary, we recommend performing an 11-point dose curve for each compound using at least *duplicate* wells for each dilution. The protocol and volumes described below are designed for a complete 96-well plate.



## DAY 2 OR 3: ANTAGONIST COMPOUND PREPARATION AND ADDITION

1. Dissolve antagonist compound in the vehicle of choice (DMSO, Ethanol, PBS or other) at the desired concentration.
2. Prepare 3-fold serial dilutions of antagonist compound in CP Reagent containing the appropriate solvent (DMSO, ethanol, PBS or other). The concentration of each dilution should be prepared at **22X** of the final screening concentration (i.e. 5 µL antagonist compound will be used in a final volume of 110 µL). For each dilution, the final concentration of solvent should remain constant.

### Preparation of 11-point dose curve serial dilutions:

We recommend starting with a concentration that is **50X** the expected  $IC_{50}$  value for the compound (**1100X**  $IC_{50}$  would be the final screening concentration).

**Example:** If the expected  $IC_{50}$  is 10 nM, prepare the highest working concentration at 11 µM.