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
Protocol for the binding assay of SARS-CoV-2 Spike glycoprotein and ACE2

Code No. HAK-SPD_bio-1

Code No. HAK-SPD_UL-1

Code No. HAK-ACE2_UL-1

August 1, 2020

【Experiment 1, binding assay】

Reagents

- Human ACE2, His Tag
- Biotin labeled SARS-CoV-2 Spike glycoprotein, Rabbit IgG1 Fc tag
- HRP-labeled streptavidin
- Immobilization buffer : PBS (-)
- Dilution buffer/Blocking buffer : 1%BSA/ PBS (-)
- Washing buffer : PBS(-) containing 0.05%Tween20
- HRP substrate : TMB etc.
- Stopping solution : 2N H₂SO₄


- 96-well plate
- Plate seal
- Plate shaker
- Plate reader

Preparation of reagents

1. Dilute human ACE2, His Tag with PBS(-) to 1 μg /mL. 100 μL of the solution is needed per well.
2. Dilute biotin labeled SARS-CoV-2 Spike glycoprotein, Rabbit IgG1 Fc tag with 1%BSA/PBS(-) to 29-21000ng/mL (serial dilution by three-fold). 100 μL of the solution is needed per well.
3. Dilute HRP-labeled streptavidin according to manufacturer's manual. 100 μL of the solution is needed per well.

Method for the measurement

1. Dilute human ACE2, His Tag.
2. Add 100 μL /well of human ACE2 in 96-well plate.
3. Stand at 4°C overnight.
4. Remove the reagent completely and add 300 μL /well of washing buffer. Discard the washing buffer and repeat this washing process twice.
5. Add 200 μL /well of 1%BSA/PBS(-).
6. As a blocking procedure, stand the plate at room temperature for 2 hours or at 4°C overnight.
7. Remove the reagent completely and add 300 μL /well of washing buffer. Discard the washing buffer and repeat this washing process twice.
8. Add 100 μL /well of diluted biotin labeled SARS-CoV-2 Spike glycoprotein, Rabbit IgG1 Fc tag.
9. Stand at room temperature for 2 hours.
10. Remove the reagent completely and add 300 μL /well of washing buffer. Discard the washing buffer and repeat this washing process twice.
11. Add 100 μL /well of diluted HRP-labeled streptavidin.
12. Stand at room temperature (reaction time should be referred to the manufacturer's manual) .
13. Remove the reagent completely and add 300 μL /well of washing buffer. Discard the washing buffer and repeat this washing process twice.
14. Add 100 μL /well of substrate solution to react at room temperature.
15. When the color developed, add 50 μL /well of stopping solution.
16. Measure the optical density (OD450) of each well by plate reader.

【Experiment 2 : Screening for the binding inhibitors by competition assay】

Reagents

- Human ACE2, His Tag
 - Bioin labeled SARS-CoV-2 Spike glycoprotein, Rabbit IgG1 Fc tag
 - SARS-CoV-2 Spike glycoprotein, Rabbit IgG1 Fc tag
 - HRP-labeled streptavidin
 - Immobilization buffer : PBS (-)
 - Dilution buffer/Blocking buffer : 1%BSA/ PBS (-)
 - Washing buffer : PBS(-) containing 0.05%Tween20
 - HRP substrate : TMB etc.
 - Stopping solution : 2N H₂SO₄
-
- 96-well plate
 - Plate seal
 - Plate shaker
 - Plate reader


Preparation of reagents

1. Dilute human ACE2, His Tag with PBS(-) to 1μg /mL. 100 μL of the solution is needed per well.
2. Dilute biotin labeled SARS-CoV-2 Spike glycoprotein, Rabbit IgG1 Fc tag with 1%BSA/PBS(-) to 1200ng/mL. 50 μL of the solution is needed per well to make final concentration of 600ng/mL.
3. Dilute SARS-CoV-2 Spike glycoprotein, Rabbit IgG1 Fc tag with 1%BSA/PBS(-) to 20-10000ng/mL (serial dilution by two-fold) to make final concentration of 10-5000ng/mL.
4. Dilute HRP-labeled streptavidin according to manufacturer's manual.

Method for the measurement

1. Dilute human ACE2, His Tag.
2. Add 100 μL/well of human ACE2 in 96-well plate.
3. Stand at 4°C overnight.
4. Remove the reagent completely and add 300μL/well of washing buffer. Discard the washing buffer and repeat this washing process twice.
5. Add 200μL /well of 1%BSA/PBS(-).
6. As a blocking procedure, stand the plate at room temperature for 2 hours or at 4°C

- overnight.
7. Remove the reagent completely and add 300 μ L/well of washing buffer. Discard the washing buffer and repeat this washing process twice.
 8. Add 50 μ L/well of diluted biotin labeled SARS-CoV-2 Spike glycoprotein, Rabbit IgG1 Fc tag.
 9. Add 50 μ L/well of SARS-CoV-2 Spike glycoprotein, Rabbit IgG1 Fc tag as a control inhibitor or any inhibitor samples to be tested.
 10. Stand at room temperature for 2 hours.
 11. Remove the reagent completely and add 300 μ L/well of washing buffer. Discard the washing buffer and repeat this washing process twice.
 12. Add 100 μ L/well of diluted HRP-labeled streptavidin.
 13. Stand at room temperature (reaction time should be referred to the manufacturer's manual).
 14. Remove the reagent completely and add 300 μ L/well of washing buffer. Discard the washing buffer and repeat this washing process twice.
 15. Add 100 μ L/well of substrate solution to react at room temperature.
 16. When the color developed, add 50 μ L/well of stopping solution.
 17. Measure the optical density (OD450) of each well by plate reader.
 18. Check for the decrease in the absorbance caused by the binding inhibitors.

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