

Ghent, March 12, 2018

CUSTOMER INFORMATION LETTER

Subject:

INNO-BIA AlzBio3 <96T, RUO>; product code 80584

Dear Valued Customer,

We hereby would like to inform you that as of the next batch (batch number 405221), the composition of the INNO-BIA AlzBio3 assay kit will change, i.e. the Coated Beads 100x vial, which contained a mix of TAU, A β (1-42) and pTAU(181P) beads, will be replaced by three separate vials containing each one type of coated beads (TAU, A β (1-42) and pTAU(181P) respectively). Please check the Package Insert for further details.

This change was done to improve the performance of the INNO-BIA AlzBio3 product by eliminating the possibility of non-specific reactions and hence increasing product stability.

Please do not hesitate to contact your local Account Manager or our Customer Support department (customer.support@fujirebio.com; tel. +32 9 329 16 11; fax. +32 9 329 17 75) in case you have any further questions.

Yours sincerely,

Nathalie Le Bastard
EMEA Marketing Manager Neuro

