

マルチゲル® II 分析範囲

商品名	品番	ウェル数	ゲル濃度	分析範囲/ SDS-PAGE (kDa)																	
				3	4	10	12	15	20	30	35	45	85	100	130	150	200	250	450	500	
グラフィエント	Multi Gel II mini 2/15	414855	13	2~15%	30~500 kDa																
		414862	17																		
	Multi Gel II mini 4/20	414879	13	4~20%	15~250 kDa																
		414886	17																		
	Multi Gel II mini 5/10	441776	13	5~10%	35~450 kDa																
		443114	17																		
	Multi Gel II mini 8/16	417269	13	8~16%	20~200 kDa																
		417276	17																		
Multi Gel II mini 10/20	414893	13	10~20%	12~130 kDa																	
	414909	17																			
Multi Gel II mini 2D-10/20	415074	1	10~20%	12~130 kDa																	
Multi Gel II mini 15/20	432026	13	15~20%	3~85 kDa																	
	443121	17																			
Multi Gel II mini 15/25	414916	13	15~25%	3~85 kDa																	
	414923	17																			
均一ゲル	Multi Gel II mini 5	443138	13	5%	100~500 kDa																
		443145	17																		
	Multi Gel II mini 7.5	414930	13	7.5%	45~250 kDa																
		414947	17																		
	Multi Gel II mini 10	414954	13	10%	30~200 kDa																
414961		17																			
Multi Gel II mini 12.5	414978	13	12.5%	20~150 kDa																	
	414985	17																			
Multi Gel II mini 15	443152	13	15%	10~150 kDa																	
	443169	17																			
ナローレンジ	Multi Gel II mini 6/9	414992	13	6~9%	45~250 kDa																
		415005	17																		
	Multi Gel II mini 9/11	415012	13	9~11%	30~200 kDa																
		415029	17																		
Multi Gel II mini 11/14	415036	13	11~14%	20~150 kDa																	
	415043	17																			
Multi Gel II mini 14/16	415050	13	14~16%	15~100 kDa																	
	415067	17																			