



Hooke Laboratories

Tuesday, April 13 2021

Dear customers,

The syringes we normally use for our kits are now being used to deliver COVID vaccines. Unfortunately, across the US, the supply of these syringes has been exhausted. During this shortage we will be using different tuberculin syringes.

The new syringes have an advantage that we can place more of them into a box, which in many cases will result in lower shipping cost. Once the syringes that we previously used are available again, we will consider switching back.

Rarely, we have also found that MOG₃₅₋₅₅ emulsion can block needles. Moving forward, we have switched to 23-gauge needles in all MOG₃₅₋₅₅ emulsion kits, which helps alleviate blocked needles in these instances.

We would appreciate your feedback on these new syringes and needles, while we continue internal testing.

Thank you for your support and understanding.

Best regards,

Kacey Caradonna, Ph.D.
Director of Product Development
k.caradonna@hookelabs.com