

Manufactured By:


Leinco Technologies, Inc.
410 Axminister Dr.
St. Louis, Missouri 63026 USA
Phone: 636-230-9477
Fax: 636-527-5545

Product Information/Name: Anti-Mouse PD-L1 Purified Low Endotoxin Functional Grade, Monoclonal Antibody, PLATINUM

Product Code: P371
Lot Number: See Lot Specific Datasheet
Parent Lot: See Lot Specific Datasheet
Clone: 10F.9G2
Isotype: Rat IgG_{2b} κ
Storage: -2°C to -8°C

Formulation: This monoclonal antibody is aseptically packaged and formulated in 0.01 M phosphate buffered saline (PBS) pH 7.2-7.4, 150 mM NaCl with no carrier protein, potassium or preservatives added.

Test	Specification	Result
Concentration	≥ 1.0 mg/ml	Conforms
Endotoxin	<0.5 EU/mg	Conforms
Purity by SDS-PAGE	≥95%	Conforms
Aggregation Level	≤1%	Conforms



Jennifer Kable
Quality Assurance

2/14/2018

Date: