



dvmedia Reference Sheet



1 Fibroblast Cellutions™Media			Appropriate Cell Types		Cell Sku
Media Sku	I-Gro-001	Sizes Available	100 mL	A Skin Fibroblasts (pre natal)	pI001-f
			500 mL	B Skin Fibroblasts (post natal)	aI001-f
				C Diseased Fibroblasts (ALS)	aI001-f-ALS
				D Diseased Fibroblasts (DMD)	aI001-f-DMD
				E Diseased Fibroblasts (TM)	aI001-f-TM
				F Diseased Fibroblasts (AC)	aI001-f-AC
				G Diseased Fibroblasts (GM)	aI001-f-GM
				H Diseased Fibroblasts (NF)	aI001-f-NF
				I Diseased Fibroblasts (SLE)	aI001-f-SLE
				J Diseased Fibroblasts (Ps)	aI001-f-Ps
				K Diseased Fibroblasts (GBS)	aI001-f-GBS
				L Diseased Fibroblasts (AVM)	aI001-f-AVM
				M Diseased Fibroblasts (MPS)	aI001-f-MPS
				N Diseased Fibroblasts (DT1)	aI001-f-DT1
				O Diseased Fibroblasts (RA)	aI001-f-RA
				P Diseased Fibroblasts (LCP)	aI001-f-LCP
				2 Stromal Cellutions™Media	
Media Sku	H-Gro-005	Sizes Available	100 mL	A Bone Marrow Stromal Cells (pre natal)	pH005-f
			500 mL	B Bone Marrow Stromal Cells (post natal)	aH005-f
				C Diseased Stromal Cells (ALS)	aH005-f-ALS
				D Diseased Stromal Cells (DMD)	aH005-f-DMD
				E Diseased Stromal Cells (DT2)	aH005-f-DT2
				F Bone Marrow Mononuclear Cells	aH002-f-2.5, -10 and -25
				G Diseased Bone Marrow Mononuclear Cells (ALS)	aH002-f-ALS-2.5, -10 and -25
				H Diseased Bone Marrow Mononuclear Cells (DMD)	aH002-f-DMD-2.5, -10 and -25
				I Diseased Bone Marrow Mononuclear Cells (DT2)	aH002-f-DT2-2.5, -10 and -25
				J Umbilical Cord Cells	aC006-f
				K CD 133+ Bone Marrow Cells	pH004-f
				L CD34+ Bone Marrow Cells	pH003-f
				M Kidney Cells	pU001-f
				N Tongue Cells	pD009-f
				O Stomach Cells	pD005-f
				P CD34+ Liver Cells	pD002-f
				3 Neural Cellutions™Media	
Media Sku	N-Gro-001	Sizes Available	100 mL	A Neural Cells	pN001-f
			500 mL	B Neural Progenitor Cells	pN003-f
				C PSA-NCAM+ Neural Cells	pN004-f
				D A2B5+ Neural Cells	pN006-f

4 Neural Pro-Conditioned Media™		Appropriate Cell Types	Cell Sku
Media Sku Sizes Available N-Pro-001 25 mL 50 mL 100 mL * N-Pro is designed to facilitate differentiation of stem or progenitor cells into neural cells.		A Bone Marrow Stromal Cells*	pH005-f
		B Bone Marrow Stromal Cells*	aH005-f
		C Umbilical Cord Cells *	aC006-f
		D CD 133 Positive Bone Marrow Cells*	pH004-f
		E Human CD34+ Bone Marrow Cells*	pH003-f
		F Bone Marrow Mononuclear Cells*	aH002-f-2.5, -10 and -25
		G Bone Marrow Mononuclear Cells (DT2)*	aH002-f-DT2-2.5, -10 and -25
		H Bone Marrow Mononuclear Cells (DMD)*	aH002-f-DMD-2.5, -10 and -25
		I Bone Marrow Mononuclear Cells (ALS)*	aH002-f-ALS-2.5, -10 and -25
		J Bone Marrow Stromal Cells (ALS)*	aH005-f-ALS
		K Bone Marrow Stromal Cells (DT2)*	aH005-f-DT2
		L Bone Marrow Stromal Cells (DMD)*	aH005-f-DMD
5 Epithelial Pro-Conditioned Media		Appropriate Cell Types	Cell Sku
Media Sku Sizes Available D-Pro-015 25 mL 50 mL 100 mL		A Whole Kidney Cells	pU001-f
		B Kidney Epithelial Cells	pU002-f
		C Whole Small Intestines Cells	pD007-f
		D Whole Large Intestines Cells	pD008-f
		F Small Intestines Epithelial Cells	pD015-f
		E Esophagus Epithelial Cells	pD016-f
6 Muscle Cellutions™Media		Appropriate Cell Types	Cell Sku
Media Sku Sizes Available M-Gro-001 500 mL		A Uncultured Skeletal Muscle Cells	pM001-f
		B Skeletal Muscle Progenitor Cells	pM002-f
		C Skeletal Muscle Progenitor Cells	aM002-f
		D Cultured Skeletal Muscle Cells	pM003-f
		E Cultured Skeletal Muscle Cells	aM003-f
		F Skeletal Muscle Progenitor Cells (DMD)	aM002-f-DMD
		G Human Skeletal Muscle Cells (DMD)	aM003-f-DMD
7 Muscle Cellutions™ Differentiation Media		Appropriate Cell Types	Cell Sku
Media Sku Sizes Available M-Diff-001 100 mL 500 mL		A Skeletal Muscle Progenitor Cells	pM002-f
		B Skeletal Muscle Progenitor Cells	aM002-f
		C Skeletal Muscle Progenitor Cells (DMD)	aM002-f-DMD
8 Cardiomyocyte Cellutions™ Maintenance Media		Appropriate Cell Types	Cell Sku
Media Sku Sizes Available C-Main-001 500 mL		A Cardiomyocytes	pC008-f
		B Cardiomyocytes	aC008-f
9 Cardiomyocyte Cellutions™ Differentiation Media		Appropriate Cell Types	Cell Sku
Media Sku Sizes Available C-MDiff-001 100 mL 500 mL		A Cardiac Progenitor Cells	pC015-f
		B Cardiac Progenitor Cells	aC015-f
10 Cardiac Cellutions™Media		Appropriate Cell Types	Cell Sku
Media Sku Sizes Available C-MGro-001 500 mL		A Cardiac Cells	pC001-f
		B Cardiac Stromal Cells	pC009-f
		C Cardiac Stromal Cells	aC009-f
		D Cardiac Progenitor Cells	pC015-f
		E Cardiac Progenitor Cells	aC015-f
		F Aortic Cells	pC016-f