



Technical Bulletin

Affected Product Description: UBIQAPTURE-Q® kit


Affected Product Catalog Number: BML-UW8995-0001

Replacement Product Catalog Number: BML-UW8995A-0001

Reason: New Mono- and polyubiquitinated conjugates monoclonal antibody

Standard Effective Date: November 20, 2020

Mono- and polyubiquitinated conjugates monoclonal antibody (FK2) (HRP conjugate) (BML-PW0150) has been replaced by Mono- and polyubiquitinated conjugates mAb (UBCJ2) (HRP conjugate) (ENZ-ABS840HRP).



Western Blot analysis of ubiquitin enrichment of lysate derived ubiquitinated proteins. Ubiquitin enrichment experiment set-up and run as described in UBIQAPTURE-Q® kit. Ubiquitin-protein conjugates present in starting material, unbound fraction and elution fraction were detected by western blotting using the provided ubiquitin-conjugate specific HRP-linked antibody. Left - Detection of ubiquitin-protein conjugates with BML-PW0150, provided in kit Prod. No. BML-UW8995. Right - Detection of ubiquitin-protein conjugates with ENZ-ABS840HRP, provided in kit Prod. No. BML-UW8995A. Key: SM = Starting Material, UF = Unbound Fraction, WA = Wash Fraction, and EL = Elution Fraction.