

■ Human IFN Alpha

品番	商品	詳細	精製方法	エンドトキシン レベル	有効期限
11200-1	Universal Type I Interferon	Alpha Interferon hybrid (human Interferon-Alpha A/D [BglII]) constructed from recombinant Hu-IFN-alpha A and Hu-IFN-alpha D.	Purified from E. coli by a combination of ion exchange, hydrophobic interaction and size exclusion chromatography	< 1 EU/μg	納品後1年
11200-2	Universal Type I Interferon	Alpha Interferon hybrid (human Interferon-Alpha A/D [BglII]) constructed from recombinant Hu-IFN-alpha A and Hu-IFN-alpha D.	A combination of ion exchange, hydrophobic interaction and size exclusion chromatography.	< 1 EU/μg	納品後1年
11002-1	Human Interferon Alpha Sampler Set	Recombinant Human Interferon-Alpha Sampler, 12 individual subtypes including: Human Interferon-Alpha A, B2, C, D, F, G, H2, I, J1, K, 4b, WA	Purified from E. coli by a combination of ion exchange, hydrophobic interaction and size exclusion chromatography	< 1 EU/μg	納品後1年
11100-1	Human Interferon Alpha A (Alpha 2a)	Recombinant Human Interferon-Alpha A (Alpha 2a), Gene IFNA2	A combination of ion exchange, hydrophobic interaction and size exclusion chromatography.	< 1 EU/μg	納品後1年
11101-1	Human Interferon Alpha A (Alpha 2a), carrier-free	Recombinant Human Interferon-Alpha A (Alpha 2a), Gene Interferon-A2	A combination of ion exchange, hydrophobic interaction and size exclusion chromatography.	< 1 EU/μg	納品後1年
11101-2	Human Interferon Alpha A (Alpha 2a), carrier-free	Recombinant Human Interferon-Alpha A (Alpha 2a), Gene Interferon-A2	A combination of ion exchange, hydrophobic interaction and size exclusion chromatography.	< 1 EU/μg	納品後1年
11105-1	Human Interferon Alpha 2 (Alpha 2b)	Recombinant Human Interferon-Alpha 2 (Alpha 2b), Gene Interferon-A2	A combination of ion exchange, hydrophobic interaction and size exclusion chromatography.	< 1 EU/μg	納品後1年
11115-1	Human Interferon Alpha B2 (Alpha 8)	Recombinant Human Interferon-Alpha B2 (Alpha 8), Gene Interferon-A8	Purified from E. coli by a combination of ion exchange, hydrophobic interaction and size exclusion chromatography	< 1 EU/μg	納品後1年
11120-1	Human Interferon Alpha C (Alpha 10)	Recombinant Human Interferon-Alpha C (Alpha 10), Gene Interferon-A10	Purified from E. coli by a combination of ion exchange, hydrophobic interaction and size exclusion chromatography	< 1 EU/μg	納品後1年
11125-1	Human Interferon Alpha D (Alpha 1)	Recombinant Human Interferon-Alpha D (Alpha 1 [Val 114]), Gene Interferon-A1	Purified from E. coli by a combination of ion exchange, hydrophobic interaction and size exclusion chromatography	< 1 EU/μg	納品後1年

品番	商品	詳細	精製方法	エンドトキシン レベル	有効期限
11130-1	Human Interferon Alpha F (Alpha 21)	Recombinant Human Interferon-Alpha F (Alpha 21), Gene Interferon-A21	Purified from E. coli by a combination of ion exchange, hydrophobic interaction and size exclusion chromatography	< 1 EU/μg	納品後1年
11135-1	Human Interferon Alpha G (Alpha 5)	Recombinant Human Interferon-Alpha G (Alpha 5), Gene Interferon-A5	Purified from E. coli by a combination of ion exchange, hydrophobic interaction and size exclusion chromatography	< 1 EU/μg	納品後1年
11145-1	Human Interferon Alpha H2 (Alpha 14)	Recombinant Human Interferon-Alpha H2 (Alpha 14), Gene Interferon-A14	Purified from E. coli by a combination of ion exchange, hydrophobic interaction and size exclusion chromatography	< 1 EU/μg	納品後1年
11150-1	Human Interferon Alpha I (Alpha 17)	Recombinant Human Interferon-Alpha I (Alpha 17), Gene Interferon-A17	Purified from E. coli by a combination of ion exchange, hydrophobic interaction and size exclusion chromatography	< 1 EU/μg	納品後1年
11160-1	Human Interferon Alpha J1 (Alpha 7)	Recombinant Human Interferon-Alpha J1 (Alpha 7), Gene Interferon-A7	Purified from E. coli by a combination of ion exchange, hydrophobic interaction and size exclusion chromatography	< 1 EU/μg	納品後1年
11165-1	Human Interferon Alpha K (Alpha 6)	Recombinant Human Interferon-Alpha K (Alpha 6), Gene Interferon-A6	Purified from E. coli by a combination of ion exchange, hydrophobic interaction and size exclusion chromatography	< 1 EU/μg	納品後1年
11175-1	Human Interferon Alpha 1 (Alpha D)	Recombinant Human Interferon-Alpha 1 (Alpha D [Ala 114]), Gene Interferon-A1	Purified from E. coli by a combination of ion exchange, hydrophobic interaction and size exclusion chromatography	< 1 EU/μg	納品後1年
11177-1	Human Interferon Alpha 4a (Alpha M1)	Recombinant Human Interferon-Alpha 4a (Alpha M1), Gene Interferon-A4	Purified from E. coli by a combination of ion exchange, hydrophobic interaction and size exclusion chromatography	< 1 EU/μg	納品後1年
11180-1	Human Interferon Alpha 4b (Alpha 4)	Recombinant Human Interferon-Alpha 4b (Alpha 4), Gene Interferon-A4	Purified from E. coli by a combination of ion exchange, hydrophobic interaction and size exclusion chromatography	< 1 EU/μg	納品後1年
11190-1	Human Interferon Alpha WA (Alpha 16)	Recombinant Human Interferon-Alpha WA (Alpha 16), Gene Interferon-A16	Purified from E. coli by a combination of ion exchange, hydrophobic interaction and size exclusion chromatography	< 1 EU/μg	納品後1年