

**Manufactured By:**

Leinco Technologies, Inc.  
410 Axminister Dr.  
St. Louis, Missouri 63026 USA  
Phone: 636-230-9477  
Fax: 636-527-5545

**Product Information/Name:** Anti-Mouse IFNAR1 Purified Low Endotoxin Functional Grade, Monoclonal Antibody, PLATINUM

Product Code: I-1188  
Lot Number: See Lot Specific Datasheet  
Parent Lot: See Lot Specific Datasheet  
Clone: MAR1-5A3  
Isotype: Mouse IgG<sub>1</sub>  
Storage: -2°C to -8°C

**Formulation:** This monoclonal antibody is aseptically packaged and formulated in 0.01 M phosphate buffered saline (PBS) pH 7.2-7.4, 150 mM NaCl with no carrier protein, potassium or preservatives added.

Test	Specification	Result
Concentration	≥ 1.0 mg/ml	Conforms
Endotoxin	<0.5 EU/mg	Conforms
Purity by SDS-PAGE	≥95%	Conforms
Aggregation Level	≤1%	Conforms

  
Jennifer Kable  
Quality Assurance

  
Date: