

■ 哺乳類 IFNs（マウス/ヒト以外）

品番	商品	詳細	精製方法	エンドトキシン レベル	有効期限
16105-1	Cynomolgus Interferon Alpha 2 (Ile 16), mammalian	Recombinant Cynomolgus Interferon-Alpha 2 (Ile 16), mammalian expressed	Purified from mammalian tissue culture media by a combination of ion exchange, hydrophobic interaction and size exclusion chromatography.	< 1 EU/μg	製造から1年
14110-1	Rhesus/Cynomolgus Interferon Alpha 2	Recombinant Rhesus Interferon-Alpha 2		< 1 EU/μg	納品後1年
15100-1	Feline Interferon Alpha	Recombinant Cat Interferon-Alpha	A combination of ion exchange, hydrophobic interaction and size exclusion chromatography.	< 1 EU/μg	納品後1年
17105-1	Porcine Interferon Alpha	Recombinant Pig Interferon-Alpha	Purified from mammalian tissue culture media by a combination of ion exchange, hydrophobic interaction and size exclusion chromatography.	< 1 EU/μg	納品後1年
19616-1	Bovine Interferon Tau 2, mammalian	Recombinant Bovine Interferon-Tau expressed in mammalian cells. Active on cow, sheep, human and mouse cells.		< 1 EU/μg	納品後1年