

■ Mouse IFN Alpha

品番	商品	詳細	精製方法	エンドトキシン レベル	有効期限
12100-1	Mouse Interferon Alpha A	Recombinant Mouse Interferon-Alpha A	A combination of ion exchange, hydrophobic interaction and size exclusion chromatography.	< 1 EU/μg	納品後1年
12105-1	Mouse Interferon Alpha 1	Recombinant Mouse Interferon-Alpha 1	A combination of ion exchange, hydrophobic interaction and size exclusion chromatography.	< 1 EU/μg	納品後1年
12115-1	Mouse Interferon Alpha 4, mammalian	Recombinant Mouse Interferon-Alpha 4	Purified from mammalian tissue culture media by a combination of ion exchange, hydrophobic interaction and size exclusion chromatography.	< 1 EU/μg	納品後1年
12125-1	Mouse Interferon Alpha 11	Recombinant Mouse Interferon-Alpha 11	A combination of ion exchange, hydrophobic interaction and size exclusion chromatography.	< 1 EU/μg	納品後1年
12130-1	Mouse Interferon Alpha 13	Recombinant Mouse Interferon-Alpha 13	A combination of ion exchange, hydrophobic interaction and size exclusion chromatography.	< 1 EU/μg	納品後1年