

■ Mouse IFN Beta

品番	商品	詳細	精製方法	エンドキシン レベル	有効期限
12405-1	Mouse Interferon Beta, mammalian	Recombinant Mouse Interferon-Beta, mammalian expressed	Purified from mammalian tissue culture media by a combination of ion exchange, hydrophobic interaction and size exclusion chromatography.	< 1 EU/μg	製造から 24カ月
12410-1	Mouse Interferon Beta, mammalian, carrier-free	Recombinant Mouse Interferon-Beta, mammalian expressed, carrier-free	Purified from mammalian tissue culture media by a combination of ion exchange, hydrophobic interaction and size exclusion chromatography.	< 1 EU/μg	製造から 24カ月
12400-1	Mouse Interferon Beta	Recombinant Mouse Interferon-Beta	A combination of ion exchange, hydrophobic interaction and size exclusion chromatography.	< 1 EU/μg	製造から 24カ月
12401-1	Mouse Interferon Beta, carrier-free	Recombinant Mouse Interferon-Beta, carrier-free	A combination of ion exchange, hydrophobic interaction and size exclusion chromatography.	< 1 EU/μg	製造から 24カ月