

■ Mouse IFN Alpha ELISA

品番	品名	詳細	サイズ	サンプル量	測定範囲	有効期限
42115-1	VeriKine-HS Mouse Interferon Alpha All Subtype ELISA Kit	Quantitative determination of Mu-IFN Alpha in serum, plasma, and tissue culture media	1 x 96	100 µl	2.38-152 pg/ml	製造から12カ月
42120-1	VeriKine Mouse Interferon Alpha ELISA Kit	Quantitative determination of Mu-IFN Alpha in cell culture and media	1 x 96	100 µl	12.5-400 pg/ml	製造から21カ月
42120-2	VeriKine Mouse Interferon Alpha ELISA Kit	Quantitative determination of Mu-IFN Alpha in cell culture and media	5 x 96	100 µl	12.5-400 pg/ml	製造から21カ月