

■ Mouse IFN Alpha 抗体

品番	品名	詳細	抗原	種由来	アイソタイプ	免疫動物	エンドキシン レベル	純度	精製方法	バッファー	タイプ	サイズ	単位	アプリケーション	推奨されるアプリケーション ※	有効期限	クローン番号
22100-1	Anti-Mouse Interferon Alpha, Clone RMMA-1 (Mab)	Neutralizes mouse interferon alpha; binds to mouse interferon alpha	IFN alpha	Mouse	IgG1	Rat	< 1 EU/μg	> 95%	A combination of affinity, ion exchange, hydrophobic interaction and size exclusion chromatography	50mM Sodium Bicarbonate; 0.2M Sodium Chloride, pH 8.0	Monoclonal Antibody	250	μg	Neutralization; ELISA	Immunoprecipitation	納品後1年	RMMA-1
22100-3	FITC Conjugated Anti-Mouse Interferon Alpha, Clone RMMA-1 (Mab)	Binds to mouse interferon alpha	IFN alpha	Mouse	IgG1	Rat				Phosphate buffered saline (PBS) containing 1% BSA and 0.05% Kathon, pH 7.2	Monoclonal Antibody	25	μg	Flow cytometry; Immunohistochemistry	-	製造後24カ月	RMMA-1
32100-1	Anti-Mouse Interferon Alpha, Rabbit Serum (PAb)	Binds to and neutralizes mouse IFN alpha	IFN alpha	Mouse		Rabbit				Neat serum diluted in phosphate buffered saline (PBS)	Polyclonal Antibody	2x10 <sup>4</sup>	units	Neutralization	ELISA; immunoprecipitation; immunohistochemistry	納品後1年	
32120-1	Anti-Mouse Interferon Alpha, Chicken IgY, Purified (PAb)	Binds to but not neutralize mouse IFN alpha	IFN alpha	Mouse	IgY	Chicken				Phosphate buffered saline (PBS) containing 0.1% BSA	Polyclonal Antibody	500	μg	Western blot; Direct Binding ELISA	-	納品後1年	