

■ Mouse IFN Beta ELISA

品番	品名	詳細	サイズ	サンプル量	測定範囲	有効期限
42410-1	VeriKine-HS Mouse Interferon Beta Serum ELISA Kit	Quantitative determination of Mu-IFN Beta in serum, plasma and tissue culture media	1 x 96	100 µl	0.94-60 pg/ml	製造から12カ月
42410-2	VeriKine-HS Mouse Interferon Beta Serum ELISA Kit	Quantitative determination of Mu-IFN Beta in serum, plasma and tissue culture media	5 x 96	100 µl	0.94-60 pg/ml	製造から12カ月
42400-1	VeriKine Mouse Interferon Beta ELISA Kit	Quantitative determination of Mu-IFN Beta in cell culture, media and serum	1 x 96	100 µl	15.6-1000 pg/ml	製造から12カ月
42400-2	VeriKine Mouse Interferon Beta ELISA Kit	Quantitative determination of Mu-IFN Beta in cell culture, media and serum	5 x 96	100 µl	15.6-1000 pg/ml	製造から12カ月