

■ Mouse IFN Beta 抗体

品番	品名	詳細	抗原	種由来	アイソタイプ	免疫動物	エンドキシン レベル	純度	精製方法	バッファー	タイプ	サイズ	単位	アプリケーション	推奨されるアプリケーション※	有効期限	クローン番号
22400-1	Anti-Mouse Interferon Beta, Clone RMMB-1 (MAb)	Binds to mouse interferon beta; weakly neutralizes mouse interferon beta	IFN Beta	Mouse	IgG	Rat				Phosphate - buffered saline (PBS) containing 0.1% bovine serum albumin (BSA)	Monoclonal Antibody	250	μg	Direct Binding ELISA	Western blot; flow cytometry; immunoprecipitation; immunohistochemistry	納品後1年	RMMB-1
32400-1	Anti-Mouse Interferon Beta, Rabbit Serum (PAb)	Binds to and neutralizes mouse IFN beta	IFN Beta	Mouse		Rabbit				Neat serum diluted in phosphate buffered saline (PBS)	Polyclonal Antibody	2x10 <sup>4</sup>	units	Neutralization	Western blot; immunoprecipitation; immunohistochemistry; ELISA	納品後1年	
32401-1	Anti-Mouse Interferon Beta, Rabbit IgG, Protein A purified (PAb)	Binds to and neutralizes mouse IFN beta	IFN Beta	Mouse		Rabbit				0.3 M Sodium Bicarbonate; 0.2 M Sodium Chloride	Polyclonal Antibody	2x10 <sup>4</sup>	units	Neutralization	Western blot; immunoprecipitation; immunohistochemistry; ELISA	納品後1年	