

■ 3本鎖糖鎖

品名	略語	標識	メーカー	品番	包装
Asialo-, agalacto-, triantennary	NGA3	Unlabeled Unlabeled 2-AA 2-AA 2-AB 2-AB	LUD PZM LUD PZM LUD PZM	CN-NGA3-20U GKC-005300 CAA-NGA3-01 GKSA-307 CAB-NGA3-01 GKSB-307	20 μ g 20 μ g 100pM 100pM 100pM 100pM
Asialo-, agalacto-, triantennary with bisecting GlcNAc (1)	—	2-AP	MCI	PA-056	100pM
Asialo-, galactosylated triantennary	NA3	Unlabeled 2-AA 2-AB	LUD	CN-NA3-20U CAA-NA3-01 CAB-NA3-01	20 μ g 100pM 100pM
Asialo-, galactosylated triantennary (1)	NA3	Unlabeled 2-AA 2-AB 2-AP	PZM PZM PZM MCI	GKC-035300 GKSA-308 GKSB-308 PA-013	10 μ g 100pM 100pM 100pM
Asialo-, galactosylated triantennary (2)	NA3	2-AP	MCI	PA-014	100pM
Tri-sialylated-, galactosylated triantennary	A3	Unlabeled 2-AB	PZM	GKC-335300 GKSB-314	10 μ g 100pM
Tri-sialylated-, galactosylated triantennary	A3	Unlabeled 2-AA 2-AB	LUD	CN-A3-20U CAA-A3-01 CAB-A3-01	20 μ g 100pM 100pM