






■ハイブリッド型糖鎖

品名	略語	標識	メーカー	品番	包装
Hybrid-type with bisecting GlcNAc	HYBR	Unlabeled Unlabeled	LUD PZM	CN-HYBRID-20U GKH003510	20 μ g 10 μ g
	2-AA 2-AB 2-AB 2-AP	LUD LUD PZM MCI	CAA-HYBRID-01 CAB-HYBRID-01 GKSB-111 PA-055	100pM 100pM 100pM 100pM	
Standard Sugar Chain PA-001	—	2-AP	MCI	PA-001	100pM
					
Standard Sugar Chain PA-066	—	2-AP	MCI	PA-066	100pM
					
Standard Sugar Chain PA-067	—	2-AP	MCI	PA-067	100pM
					
Standard Sugar Chain PA-068	—	2-AP	MCI	PA-068	100pM
					
Standard Sugar Chain PA-541	—	2-AP	MCI	PA-541	100pM
