

■ その他動物種 IFNs 抗体

品番	品名	詳細	抗原	種由来	アイソタイプ	免疫動物	エンドキシン レベル	純度	精製方法	バッファー	タイプ	サイズ	単位	アプリケーション	推奨されるアプリケーション※	有効期限	クローン番号
27100-1	Anti-Pig Interferon Alpha, Clone K9 (MAb)	Binds to and neutralizes pig IFN alpha	IFN alpha	Pig	IgG1	Mouse	< 1 EU/μg	> 95%	Affinity chromatography via Protein G		Monoclonal Antibody	500	μg	-	Western blot; immunoprecipitation; immunohistochemistry; ELISA	納品後1年	K9
27105-1	Anti-Pig Interferon Alpha, Clone F17 (MAb)	Binds to and neutralizes pig IFN alpha	IFN alpha	Pig	IgG1	Mouse	< 1 EU/μg	> 95%	Affinity chromatography via Protein G		Monoclonal Antibody	500	μg	-	Western blot; immunoprecipitation; immunohistochemistry; ELISA	納品後1年	F17
33100-1	Anti-Rat Interferon Alpha, Rabbit IgG, Protein A purified (PAb)	Neutralizes both natural and recombinant rat interferon alpha.	IFN alpha	Rat	IgG	Rabbit			Ammonium sulfate precipitation and protein - A affinity chromatography	0.5 ml Phosphate - buffered saline (PBS) + 125 mM trehalose prior to lyophilization	Polyclonal Antibody	500	μg	-	Western blot; In vitro neutralization	-	
33400-1	Anti-Rat Interferon Beta, Rabbit IgG, Protein A purified (PAb)	Neutralizes both natural and recombinant rat interferon beta	IFN Beta	Rat		Rabbit			Antibodies were sequentially purified by a mmonium sulphate precipitation and Protein A - affinity chromatography	0.5 ml Phosphate - buffered saline (PBS) + 125 mM trehalose prior to lyophilization	Polyclonal Antibody	500	μg	-	Western blot; ELISA; Neutralization	製造後36カ月	