

■ Rat IFNs

品番	商品	詳細	精製方法	エンドトキシン レベル	有効期限
13100-1	Rat Interferon Alpha	Recombinant Rat Interferon-Alpha			製造から 36カ月
13400-1	Rat Interferon Beta	Recombinant Rat Interferon-Beta. Active on Rat and Mouse cells.			製造から 36カ月
13500-1	Rat Interferon Gamma	Recombinant Rat Interferon-Gamma	Purified from E. coli by a combination of ion exchange, hydrophobic interaction and size exclusion chromatography	< 1 EU/μg	納品後1年
13101-1	Rat Interferon Alpha 1	Recombinant Rat Interferon Alpha 1 expressed in mammalian cells			納品後1年
13114-1	Rat Interferon Alpha 14	Recombinant Rat Interferon Alpha 14 expressed in mammalian cells			納品後1年