

製品使用文献

品番：4BDX-18K8



• Oh EY, et al. Comparison of fertility competence in toll-like receptor 4 (TLR4)-knock out male mice fed a high-fat diet. *J Biomed Transl Res.* 2023; 24(2):41-52

• Boersma, et al. Influence of sperm cryopreservation on sperm motility and proAKAP4 concentration in mice. *Reprod Med Biol.* Vol21,e12480. 2022

• Delehedde M, Carracedo S, Duchene B, Demouveaux B, Remy G, Kervoaze G, Aubry S, Gevaert MH, Maurage CA, Jonckheere N, Van Seuningen I, Philippe Gosset G, Desseyn JL, Sergeant N and Pichavant M .Cigarette smoking affects the sperm-specific proAKAP4 concentrations and impairs both spermatogenesis and sperm quality. *Andrology.* Vol8,Supp 1,45-46. 2020

• Sergeant , et al. Proteolysis of proAKAP4 in semen as a regulatory sensor of sperm quality and functionality. *Andrology.* Vol8,Supp 1,44-45. 2020

• Sergeant, et al. The sperm specific protein proAKAP4 as an innovative marker to evaluate sperm quality and fertility. *Journal of Dairy & Veterinary Sciences.* Vol11,43466. 2019

• Sergeant, et al. Investigating proteomic methods and tools to assess sperm quality. *Animal Reproduction Science.* Vol169, 99–135. 2016