

製品使用文献

品番：4VDX-18K3



- Carracedo , et al. ProAKAP4 protein marker: towards a functional approach to male fertility. *Animal Reproduction Science*. 2022 (in press)
- Dordas-Perpinya , et al. ProAKAP4 concentrations in stallion and donkey: Comparison with kinetic parameters, motility and concentration. *Reproduction in Domestic Animals*. Vol57,1,61. 2022
- Dordas-Perpinya , et al. ProAKAP4 correlation with motion parameters and spermatozoa subpopulations in stallion. *Reproduction in Domestic Animals*. 2022
- Dordas-Perpinya , et al. Characterization of proAKAP4 expression, metabolism and localization in Equus asinus spermatozoa. *Reproduction, Fertility and Development*. 2022
- Blommaert, et al. First results about ProAKAP4 concentration in stallion semen after cryopreservation in two different freezing media. *Cryobiology*. Vol102,133-35 Oct. 2021
- Sergeant , et al. Proteolysis of proAKAP4 in semen as a regulatory sensor of sperm quality and functionality. *Andrology*. Vol8,Supp 1,44-45. 2020
- Carracedo , et al. The sperm specific proAKAP4 polypeptide exhibited conserved functions, localizations and metabolism among mammals. *Animal Reproduction Science*. Vol220,106448. 2020
- Sergeant , et al. Assessment of the proAKAP4 concentrations as a pertinent sperm parameter to select extenders and preservatives of stallion semen. *Animal Reproduction Science*. Vol220,106449. 2020

- Blommaert , et al. Expression, localization, and concentration of A-kinase anchor protein 4 (AKAP4) and its precursor (proAKAP4) in equine semen: promising marker correlated to the total and progressive motility in thawed spermatozoa. *Theriogenology*. Vol131,52-60. 2019

- Sergeant, et al. The sperm specific protein proAKAP4 as an innovative marker to evaluate sperm quality and fertility. *Journal of Dairy & Veterinary Sciences*. Vol11,43466. 2019

- Blommaert, et al. The Horse 4MID® Kit as a promising tool to compare different freezing media to preserve stallion semen. *Rev. Bras. Reprod. Anim.*Vol43,2,585. 2019

- Delehedde , et al. Concentration of proAKAP4 as a pertinent read-out of sperm quality in mammals. *Animal Reproduction Science*. Vol194,24. 2018

- Blommaert , et al. Significant correlation between the proAKAP4 concentration and the total and progressive motility in stallion sperm after thawing. *Journal of Equine Veterinary Science*. Vol 6 . 6 ,43. 2018

- Sergeant, et al. Investigating proteomic methods and tools to assess sperm quality. *Animal Reproduction Science*. Vol169, 99–135. 2016

- Dordas-Perpinyà, et al. ProAKAP4 concentration is related with sperm motility and motile sperm subpopulations in frozen-thawed horse semen. *Animals*. Vol12, 3417. 2022

- Smith, et al. Preliminary Study of Effects of Centrifugation of Stallion Semen on Motility and Pro-AKAP4 Concentration over 72 Hours Post-collection. . *Journal of Agricultural Science and Technology*. Vol12. 65-72. 2022